



Aeon7200

Anaesthesia Machine

CE 0123

**AEOMED**  
Reliable Quality Thoughtful Service

- 7" color LCD screen, different parameters distinguished by colors for easy recognition
- Novel designed ventilator pneumatic system decreases the consumption of Oxygen
- Built-in NO WATER VAPOR CONDENSATION sampling device ensures ZERO maintenance
- Fast-release function by APL valve improve the efficiency of resuscitation after surgery
- Advanced one button switch of manual / auto provides convenient pneumatic / auto controlling interchanges
- Optional AGSS device for more protection
- With CE certificate, meets EU clinical safety requirement

## Technical Specification

### Physical & Environmental Specification

- Dimensions (HxLxW): 1386x839x628mm
- Weight: 85kg (basic unit)
- Storage Temperature: -20~+55°C
- Working Temperature: +10~+40°C
- Storage Humidity: ≤93%

### Electrical Specifications

- Mains: 100~240V AC, 50/60Hz
- Battery: 24V DC, minimum 120min

### Pneumatic Specifications

- Gas Supply: O<sub>2</sub>, N<sub>2</sub>O: 0.28~0.6MPa
- Flowmeter: O<sub>2</sub>: 0~1L/min; 1~10L/min  
N<sub>2</sub>O: 0~1L/min; 1~10L/min
- Gas System: Oxygen supply failure alarm  
Oxygen supply  
O<sub>2</sub> flush: 25~75L/min
- Working Mode: Closed, Semi-closed, Semi-open
- Driven Mode: Pneumatically driven and electronically controlled
- Safety Valve: ≤8kPa
- Operation mode: Man. / Vent.

### Ventilator Specification

- Patient Type: Adult, pediatric
- Ventilation Modes: IPPV, Manual;  
optional: SIPPV, SIMV, SPONT
- Setting: Touch key
- Tidal Volume: 0~1600ml
- Ventilation Frequency: 4~100bpm
- I:E Ratio: 4:1~1:8
- PEEP: 0~20 cmH<sub>2</sub>O
- Inspiratory Plateau: Off, 5%~60%Ti
- Monitoring Display: Display: 7" color wide LCD screen  
Pressure: Ppeak, Pplat, Pmean, PEEP  
Volumes: MV, Vt  
Breath rate  
Gas: FiO<sub>2</sub>  
Graph display: P-t, F-t waveforms
- Alarm: High Airway Pressure/Low Airway Pressure  
High Minute Volume/Low Minute Volume  
High FiO<sub>2</sub>/Low FiO<sub>2</sub>  
Continuous high pressure  
Power failure  
O<sub>2</sub> Supply failure  
Apnea  
Alarm Silence: ≤120 seconds

### Vaporizer

- Support 2 vaporizer (Selectatec® with interlock, optional standby vaporizer parking holder)
- Agent Type: Halothane, Enflurane, Isoflurane,  
Sevoflurane available

Remark: Above configurations include standard and option. Please check price with your Aeonmed sales representative.



**AEOMED**  
Reliable Quality Thoughtful Service

HQ: Building 9, No.26 Outer Ring West Road, Fengtai District, Beijing 100070, China  
Science Park: No. 10, Chaobai Street, Yanjiao Development Zone, Sanhe City, Hebei Province  
065201, China  
TEL: +86-10-5841 1198 FAX: +86-10-6371 8989  
Http: //www.aeonmed.com E-mail: marketing@aeonmed.com

The information contained in this leaflet is correct at the date of publication. The policy of Beijing Aeonmed Co., Ltd. is one of the continued improvement to its products. Because of this policy, Aeonmed reserves the right to make any change, which may affect the information in this leaflet without giving prior notices.

AAM7200-1902